QuantScopics® Urinalysis Microscopics Control Level 1 (Low) / Level 2 (High)

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English

Intended use

The QuantScopics Urinalysis Microscopics Control is intended for use as quality control material for the microscopic evaluation of various urine sediment methods.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Microscopic QC controls should be run each day the test is performed.¹

Summary and Explanation

Standardized microscopic evaluation of urine sediment is an important part of routine analysis of urine. Along with physical and chemical analysis, microscopic examination of urine can provide valuable information regarding not only renal and urinary tract disease, but also metabolic diseases unrelated to the kidney.² Urinary sediment microscopy generally includes the detection and identification of red blood cells, leukocytes, epithelial cells, bacteria, casts and crystals.^{3,4}

Product Description

The QuantScopics Urinalysis Microscopics Control is supplied in two levels, 4 x 120 mL bottles per level per box. It is liquid, ready-to-use, requiring no reconstitution or dilution. The control is made from human urine to which stabilized human red blood cells, human white blood cells, and uric acid, calcium oxalate and calcium phosphate crystals have been added. Preservatives have been added to maintain product integrity.

Caution

Contains human urine and human blood cells. All blood donor units comprising the human cell source material used in the manufacture of this product have been tested and found nonreactive for Hepatitis B Surface Antigen and Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.
H317 – May cause an allergic skin reaction.
P261 – Avoid breathing vapors, mist, or spray.
P280 – Wear protective gloves, protective clothing, and eye protection.
P302+P352 – IF ON SKIN: Wash with plenty of water.
P321 – Specific Treatment (see section 4 on the SDS).
P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 – Take off contaminated clothing and wash it before reuse.
P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations. Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Storage and Stability

The QuantScopics Urinalysis Microscopics Control should be stored tightly capped at 2°C–8°C when not in use. Do not freeze. When stored at 2°C–8°C the controls are stable until the expiration date stated on the label. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure

Remove the control from the refrigerator and replace the cap on the control bottle with the pour seal cap included in the control box. Allow the control to come to room temperature (18°C–25°C) for approximately 15–90 minutes depending on the volume remaining in the bottle. Mix the control thoroughly by inverting the bottle at least 20 times to assure homogeneity of the contents. Avoid foaming. Thorough mixing with each use is important in order to obtain reproducible results. Treat the control as you would a patient sample in accordance with the manufacturer's instructions for the standardized microscopic urinalysis system you are using. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the use of standardized systems in order to yield standardized, reproducible results and to enable the reporting of abnormal sediment elements per unit volume.¹ Immediately close the pour seal cap and store the control at 2°C–8°C when not in use.

Expected Ranges

The expected ranges for each type of formed element found in this lot of QuantScopics Control is listed in the table on page 4. These ranges were determined by assay of multiple bottles of the indicated lot by the methods listed. A 12 mL sample volume was used for each testing method. The samples were centrifuged at 400 RCF (relative centrifugal force) for 5 minutes. After centrifugation, urine sediment was resuspended in either ~0.5 or ~1.0 mL of remaining supernatant according to the respective cell count method procedure. The ranges listed are based on the range of elements observed in 10 high power fields. Use of other systems or protocols may yield differing results. Each laboratory should establish its own precision parameters.

Limitations

Any future changes made by the manufacturer of a test system may give different values from the indicated range. Detailed information on the limitations of each test system is included in the limitations section of the system manufacturer's package insert. Technical updates can be found on our website. The Quality Control Log can be downloaded from the Quantimetrix website at quantimetrix.com or contact Tech Support at (310) 536-0006, option 3.

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- 1 Routine Urinalysis and Collection, Transportation, and Preservation of Urine Specimens, NCCLS Documentation GP 16-A, Approved Guideline, (1995), NCCLS, Wayne PA, 19087.
- 2 Bradley GM, Benson ES, Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th ed., Philadelphia, PA: Saunders, 1974.
- 3 Bologna, CV, Understanding Laboratory Medicine, St. Louis, MO: CV Mosby Co., 1971.
- 4 Dudas HC, Lab Med, 12:765, 1981.

Microscopic Evaluation of Urine Sediment • Mikroskopische Beurteilung von Urinsediment • Evaluation microscopique du sédiment urinaire Valutazione microscopica del sedimento urinario • Evaluación microscópica de la sedimentación presente en la orina

Analytes/Method	Level 1 48301	Level 2 48302
Red Blood Cells		
AVE-764B Automated Urine Formed Elements Analyzer	0 - 19/µL	45 - 118/µL
Cen-Slide® System	0 - 12/hpf	2 - 35/hpf
FisherBrand UriSystem DeciSlide	0 - 6/hpf	0 - 20/hpf
IRIS® Diagnostics iQ200 Series Analyzers	1 - 16/µL	35 - 98/µL
Kova® Classic™ Slide 10 with grids	1 - 13/µL	16 - 89/µL
Non-grid slides (~0.5 mL resuspended sediment)	0 - 6/hpf	1 - 29/hpf
Non-grid slides (~1.0 mL resuspended sediment)	0 - 6/hpf	2 - 29/hpf
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	0 - 50/hpf	0 - 200/hpf
Slide & Coverslip (~0.5 mL resuspended sediment)	0 - 6/hpf	2 - 22/hpf
Slide & Coverslip (~1.0 mL resuspended sediment)	0 - 6/hpf	1 - 20/hpf
Sysmex UF-1000i	4 - 15/µL	58 - 84/µL
UriSed Variants / COBIO Variants	0 - 30/µL	20 - 240/µL
White Blood Cells		
AVE-764B Automated Urine Formed Elements Analyzer	0 - 10/µL	20 - 70/µL
Cen-Slide® System	0 - 11/hpf	2 - 25/hpf
FisherBrand UriSystem DeciSlide	0 - 4/hpf	0 - 18/hpf
IRIS® Diagnostics iQ200 Series Analyzers	2 - 14/µL	38 - 91/µL
Kova® Classic Slide 10 with grids	0 - 9/µL	6 - 75/µL
Non-grid slides (~0.5 mL resuspended sediment)	0 - 10/hpf	0 - 31/hpf
Non-grid slides (~1.0 mL resuspended sediment)	0 - 5/hpf	0 - 19/hpf
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	0 - 50/hpf	0 - 200/hpf
Slide & Coverslip (~0.5 mL resuspended sediment)	0 - 5/hpf	0 - 16/hpf
Slide & Coverslip (~1.0 mL resuspended sediment)	0 - 6/hpf	0 - 14/hpf
Sysmex UF-1000i	5 - 11/µL	44 - 63/µL
UriSed Variants / COBIO Variants	0 - 20/µL	20 - 100/µL
Casts		
AVE-764B Automated Urine Formed Elements Analyzer	none	none
Cen-Slide® System	none	none
IRIS® Diagnostics iQ200 Series Analyzers	none	none
Kova® Classic Slide 10 with grids	none	none
Non-grid slides (~0.5 mL resuspended sediment)	none	none
Non-grid slides (~1.0 mL resuspended sediment)	none	none
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	none	none
Slide & Coverslip (~0.5 mL resuspended sediment)	none	none
Slide & Coverslip (~1.0 mL resuspended sediment)	none	none
Sysmex UF-1000i	none	none
UriSed Variants / COBIO Variants	none	none

Analytes/Method	Level 1 48301	Level 2 48302
Crystals		
AVE-764B Automated Urine Formed Elements Analyzer	may be present	may be present
Cen-Slide® System	may be present	may be present
IRIS® Diagnostics iQ200 Series Analyzers	present	present
Kova® Classic Slide 10 with grids	may be present	may be present
Non-grid slides (~0.5 mL resuspended sediment)	may be present	may be present
Non-grid slides (~1.0 mL resuspended sediment)	may be present	may be present
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	none - present	none - present
Slide & Coverslip (~0.5 mL resuspended sediment)	may be present	may be present
Slide & Coverslip (~1.0 mL resuspended sediment)	may be present	may be present
Sysmex UF-1000i	may be present	may be present
UriSed Variants / COBIO Variants	may be present	may be present
Bacteria		
AVE-764B Automated Urine Formed Elements Analyzer	may be present	may be present
Cen-Slide® System	may be present	may be present
IRIS® Diagnostics iQ200 Series Analyzers	may be present	may be present
Kova® Classic Slide 10 with grids	may be present	may be present
Non-grid slides (~0.5 mL resuspended sediment)	may be present	may be present
Non-grid slides (~1.0 mL resuspended sediment)	may be present	may be present
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	n/a	n/a
Slide & Coverslip (~0.5 mL resuspended sediment)	may be present	may be present
Slide & Coverslip (~1.0 mL resuspended sediment)	may be present	may be present
Sysmex UF-1000i	may be present	may be present
UriSed Variants / COBIO Variants	may be present	may be present

 Quantimetrix Corporation
 2005 Manhattan Beach Blvd.
 Redondo Beach, CA 90278
 +1.310.536.0006

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This page intentionally left blank.

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MDSS GmbH