
# Dropper® A1c Diabetes Control / Level 1 & 2

## English

### Intended Use

The Quantimetrix Dropper A1c Diabetes Control is intended for the quality control of analytical procedures used to quantitate HbA1c.

### Product Description

The Quantimetrix Dropper A1c Diabetes Control is supplied in two levels as a frozen, ready-to-use material that requires no reconstitution. The control consists of human whole blood modified to target levels with reagent chemicals to produce the desired reaction when tested by the methods indicated in the Intended Use section. Stabilizers and preservatives have been added to maintain shelf life and inhibit microbial growth. Included in each kit is a lot specific control barcode card for Siemens DCA 2000/2000+Vantage users.

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units used in this product have been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods.

No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV viruses. These samples should be handled according to the Centers for Disease Control's Biosafety Level 2 recommendations.

### Storage and Stability

The Quantimetrix Dropper A1c Diabetes Control should be stored frozen at -10°C to -30°C. The control may be received thawed and should be re-frozen upon arrival. When stored frozen at -10°C to -30°C, the control remains stable until the expiration date printed on the label. Once the control is thawed and opened, it may be stored refrigerated at 2-8°C for 6 months or at room temperature at 18°C to 25°C for 21 days. Discard the control if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens according to local guidelines.

### Procedure

If frozen, allow the control to thaw for at least 30 minutes at room temperature with the cap on. Mix gently by inversion a minimum of ten times to ensure the homogeneity of the contents. Avoid foaming. Remove cap and squeeze the proper amount of control into a sample cup or on a clean flat surface for assay. Treat the control as you would a patient sample in accordance with the manufacturer's requirements for the test method. Immediately recap and return the control to its desired storage condition. Maintain the same storage condition for the life of the product. Alternating the control between refrigerated temperature and room temperature storage could adversely affect the product's shelf life. If the control is received thawed, immediately place in the desired closed vial storage condition. For further information see the storage and stability section of this insert.

### Scanning Procedure for Control Barcode Card

Siemens DCA 2000/2000+Vantage users must scan the included lot specific control barcode card before a control can be tested. Hold the card so that the barcode specific to the level of control to be tested is facing to the right. Place the card gently against the right side of the track and smoothly slide the card down. A beep sound signals a successful scan after which the corresponding control can be tested. If no beep sounds, repeat the scanning procedure.

### Expected Values

Expected values for the listed lots of controls have been established from inter-laboratory data using instrument manufacturers' reagents. Individual laboratory means should fall within the ranges listed. These values should be used as a guide in evaluating the performance of test methods. Each laboratory should establish its own ranges and specifications for the methods used to measure HbA1c. A quality control program, Quantrol®, is available from Quantimetrix website. Please inquire.

### Limitations

The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Limitations of the test method are included in the package insert for the reagent or the manual for the instrument being used. Technical updates can be found on our website.

**Note:** For all users obtaining instrument error messages for the A1c analyte, dilute the control(s) by 50% (1 part control to 1 part normal saline). Run sample as normal and record values as reported by the instrument; do not correct for the dilution.

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Quantimetrix Corporation
MDSS GmbH